

wisonic
Your needs we care



Piloter Vet Color Doppler Ultrasound A Smart Tablet Ultrasound System

Shenzhen Wisonic Medical Technology Co., Ltd.

Address: 5th Floor, No.6 Building, Pingshan Technology Park,
Taoyuan Street, Nanshan District, Shenzhen, China

Tel:+86-755-86007788

Fax:+86-755-86007799

Email: info@wisonic.cn

www.wisonic.cn



7400-00-40250

wisonic
Your needs we care

Genius Design

Piloter vet has an extremely light weight which making it more easy to move anytime and anywhere in a briefcase. With smart design and built-in customized workflow, Piloter vet provides more efficient and convenient experience for all applications.


13.3-inch touch screen which allows for horizontal or vertical viewing display



Weighing only 2Kg Extremely lightweight and portable



Complete with a built in tray and dedicated probe holders



Built in multi-port adaptor



Desk-top stand




Smart, easy to manoeuvre base



Prominent Imaging


Built on Holo platform – 64 times beamforming image processing, able to achieve up to 5000 frames/second. Piloter vet provides stunning images for all applications including animal's abdomen, reproduction as well as cardiology etc.




Cat Stomach



Dog Spleen



Dog Cardiac Color Flow



Cat Cardiac PW



Cat Kidney



Cat Intestine

Broad Set of Probes



Smart-design button on both sides of probe, make user readily interact with system via gloves-on hands.

